

Please type a plus sign (+) inside this box → ☐

PTO/SB/05 (2/98)
Approved for use through 09/30/2000. OMB 0651-0032
Patent and Trademark Office: U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

UTILITY PATENT APPLICATION TRANSMITTAL

(Only for new nonprovisional applications under 37 C.F.R. § 1.53(b))

Attorney Docket No. P-4278

First Inventor or Application Identifier Oscar J. Llorin et al.

Title Cell Disruption Method Using Sonication

Express Mail Label No. EM368495356

APPLICATION ELEMENTS

See MPEP chapter 600 concerning utility patent application contents.

1. ☒ * Fee Transmittal Form (e.g., PTO/SB/17)
(Submit an original and a duplicate for fee processing)
2. ☒ Specification [Total Pages 19]
(preferred arrangement set forth below)
 - Descriptive title of the invention
 - Cross References to Related Applications
 - Statement Regarding Fed sponsored R & D
 - Reference to Microfiche Appendix
 - Background of the invention
 - Brief Summary of the invention
 - Brief Description of the Drawings (if filed)
 - Detailed Description
 - Claim(s)
 - Abstract of the Disclosure
3. ☒ Drawing(s) (35 U.S.C. 113) [Total Sheets 1]
4. Oath or Declaration [Total Pages]
 - a. ☒ Newly executed (original or copy)
 - b. ☐ Copy from a prior application (37 C.F.R. § 1.63(d))
(for continuation/divisional with Box 17 completed)
[Note Box 5 below]
 - i. ☐ DELETION OF INVENTOR(S)
Signed statement attached deleting inventor(s) named in the prior application, see 37 C.F.R. §§ 1.63(d)(2) and 1.33(b).
5. ☐ Incorporation By Reference (useable if Box 4b is checked)
The entire disclosure of the prior application, from which a copy of the oath or declaration is supplied under Box 4b, is considered to be part of the disclosure of the accompanying application and is hereby incorporated by reference therein.

ADDRESS TO: Assistant Commissioner for Patents
Box Patent Application
Washington, DC 20231

6. ☐ Microfiche Computer Program (Appendix)
7. Nucleotide and/or Amino Acid Sequence Submission
(if applicable, all necessary)
 - a. ☐ Computer Readable Copy
 - b. ☐ Paper Copy (identical to computer copy)
 - c. ☐ Statement verifying identity of above copies

ACCOMPANYING APPLICATION PARTS

8. ☒ Assignment Papers (cover sheet & document(s))
9. ☐ 37 C.F.R. § 3.73(b) Statement (when there is an assignee) ☐ Power of Attorney
10. ☐ English Translation Document (if applicable)
11. ☐ Information Disclosure Statement (IDS)/PTO-1449 ☐ Copies of IDS Citations
12. ☐ Preliminary Amendment
13. ☒ Return Receipt Postcard (MPEP 503)
(Should be specifically itemized)
14. ☐ * Small Entity Statement(s) ☐ Statement filed in prior application, Status still proper and desired
(PTO/SB/09-12)
15. ☐ Certified Copy of Priority Document(s)
(if foreign priority is claimed)
16. ☐ Other:

NOTE FOR ITEMS 1 & 14: IN ORDER TO BE ENTITLED TO PAY SMALL ENTITY FEES, A SMALL ENTITY STATEMENT IS REQUIRED (37 C.F.R. § 1.27), EXCEPT IF ONE FILED IN A PRIOR APPLICATION IS RELIED UPON (37 C.F.R. § 1.28).

17. If a CONTINUING APPLICATION, check appropriate box, and supply the requisite information below and in a preliminary amendment:

☐ Continuation ☐ Divisional ☐ Continuation-in-part (CIP) of prior application No: _____
Prior application information: Examiner _____ Group / Art Unit: _____

18. CORRESPONDENCE ADDRESS

☐ Customer Number or Bar Code Label

(Insert Customer No. or Attach bar code label here)

or ☒ Correspondence address below

Name	Richard J. Rodrick, Esq.				
	Becton Dickinson and Company				
Address	1 Becton Drive				
	Mail Code 089				
City	Franklin Lakes	State	NJ	Zip Code	07417
Country	USA	Telephone	(201) 847-5317	Fax	(201) 84809228

Name (Print/Type)	David W. Highet, Esq.	Registration No. (Attorney/Agent)	30,265
Signature	<i>David W. Highet</i>	Date	8/3/98

Burden Hour Statement: This form is estimated to take 0.2 hours to complete. Time will vary depending upon the needs of the individual case. Any comments on the amount of time you are required to complete this form should be sent to the Chief Information Officer, Patent and Trademark Office, Washington, DC 20231. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Assistant Commissioner for Patents, Box Patent Application, Washington, DC 20231.

EXPRESS MAIL LABEL NO. EM368495356

FEE TRANSMITTAL

Patent fees are subject to annual revision on October 1.
These are the fees effective October 1, 1997.
Small Entity payments must be supported by a small entity statement,
otherwise large entity fees must be paid. See Forms PTO/SB/09-12.
See 37 C.F.R. §§ 1.27 and 1.28.

TOTAL AMOUNT OF PAYMENT (\$) 830.00

Complete if Known

Application Number	Not Yet Assigned
Filing Date	August 3, 1998
First Named Inventor	Oscar J. Llorin et al.
Examiner Name	Not Yet Assigned
Group / Art Unit	Not Yet Assigned
Attorney Docket No.	P-4278

METHOD OF PAYMENT (check one)

1. ☒ The Commissioner is hereby authorized to charge indicated fees and credit any over payments to:

Deposit Account Number 02-1666
Deposit Account Name Becton Dickinson & Co.

☐ Charge Any Additional Fee Required Under 37 C.F.R. §§ 1.16 and 1.17
☐ Charge the Issue Fee Set in 37 C.F.R. § 1.18 at the Mailing of the Notice of Allowance

2. ☐ Payment Enclosed:
☐ Check ☐ Money Order ☐ Other

FEE CALCULATION

1. BASIC FILING FEE

Large Entity Fee Code (\$)	Small Entity Fee Code (\$)	Fee Description	Fee Paid
101 790	201 395	Utility filing fee	790
106 330	206 165	Design filing fee	
107 540	207 270	Plant filing fee	
108 790	208 395	Reissue filing fee	
114 150	214 75	Provisional filing fee	
SUBTOTAL (1)			790

2. EXTRA CLAIM FEES

Total Claims 13 - 20** = 13 X 0 = 0
Independent Claims 2 - 3** = 2 X 0 = 0
Multiple Dependent 1 = 0

**or number previously paid, if greater; For Reissues, see below

Large Entity Fee Code (\$)	Small Entity Fee Code (\$)	Fee Description
103 22	203 11	Claims in excess of 20
102 82	202 41	Independent claims in excess of 3
104 270	204 135	Multiple dependent claim, if not paid
109 82	209 41	** Reissue independent claims over original patent
110 22	210 11	** Reissue claims in excess of 20 and over original patent

SUBTOTAL (2) (\$) 0

FEE CALCULATION (continued)

3. ADDITIONAL FEES

Large Entity Fee Code (\$)	Small Entity Fee Code (\$)	Fee Description	Fee Paid
105 130	205 65	Surcharge - late filing fee or oath	
127 50	227 25	Surcharge - late provisional filing fee or cover sheet	
139 130	139 130	Non-English specification	
147 2,520	147 2,520	For filing a request for reexamination	
112 920*	112 920*	Requesting publication of SIR prior to Examiner action	
113 1,840*	113 1,840*	Requesting publication of SIR after Examiner action	
115 110	215 55	Extension for reply within first month	
116 400	216 200	Extension for reply within second month	
117 950	217 475	Extension for reply within third month	
118 1,510	218 755	Extension for reply within fourth month	
128 2,060	228 1,030	Extension for reply within fifth month	
119 310	219 155	Notice of Appeal	
120 310	220 155	Filing a brief in support of an appeal	
121 270	221 135	Request for oral hearing	
138 1,510	138 1,510	Petition to institute a public use proceeding	
140 110	240 55	Petition to revive - unavoidable	
141 1,320	241 660	Petition to revive - unintentional	
142 1,320	242 660	Utility issue fee (or reissue)	
143 450	243 225	Design issue fee	
144 670	244 335	Plant issue fee	
122 130	122 130	Petitions to the Commissioner	
123 50	123 50	Petitions related to provisional applications	
126 240	126 240	Submission of Information Disclosure Stmt	
581 40	581 40	Recording each patent assignment per property (times number of properties)	40
146 790	246 395	Filing a submission after final rejection (37 CFR 1.129(a))	
149 790	249 395	For each additional invention to be examined (37 CFR 1.129(b))	
Other fee (specify) _____			
Other fee (specify) _____			

* Reduced by Basic Filing Fee Paid

SUBTOTAL (3) (\$) 40

SUBMITTED BY

Typed or Printed Name David W. Highet, Esq.

Signature David W. Highet

Date 8/2/98

Complete (if applicable)

Reg. Number 30,265

Deposit Account User ID 02-1666

Burden Hour Statement: This form is estimated to take 0.2 hours to complete. Time will vary depending upon the needs of the individual case. Any comments on the amount of time you are required to complete this form should be sent to the Chief Information Officer, Patent and Trademark Office, Washington, DC 20231. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Assistant Commissioner for Patents, Washington, DC 20231.

TITLE OF THE INVENTION

CELL DISRUPTION METHOD USING SONICATION

BACKGROUND OF THE INVENTION

5 Access to cellular components such as nucleic acids is imperative to a variety of molecular biology methodologies. Such methodologies include nucleic acid sequencing, direct detection of particular nucleic acid sequences by nucleic acid hybridization and nucleic acid sequence amplification techniques.

10 Although access to nucleic acids from the cells of some organisms does not involve particularly complex methodologies or harsh treatments, other organisms have cells from which it is particularly difficult to access nucleic acids or other cellular components. Organisms in the latter group include species of the genus *Mycobacteria*, yeast and fungi. Usually, the difficulty in cellular component access is a result of organism cell walls which are highly resistant to lysis or disruption, and/or the adherence of certain cellular components
15 such as nucleic acids to cellular proteins and other cellular substances such as pieces of cell walls.

Due to the difficulties in attempting to access nucleic acids from mycobacterial organisms the methods utilized tend to be harsh and thus not very useful with non-mycobacterial organisms. Conversely, the methods used to disrupt cells and access nucleic
20 acids from non-mycobacterial organisms are often not effective when used with mycobacterial organisms.

Two non-enzymatic methods which have been used to disrupt cells to access nucleic acids are the application of heat to cells (see U.S. Patent No. 5,376,527) and physical agitation of cells in the presence of lysogenic chemicals with or without "minibeads". For example,
25 DeWitt et al., J. Clin. Micro. 28 (11):2437-2441 (1990) describe the orbital shaking of samples containing mycobacterial cells in the presence of buffered phenol and sodium dodecyl sulfate (SDS), Hurley, S.S. et al., J. Clin. Microbiol. 25 (11) 2227-2229 (1987)

describe a combination of phenol extraction and physical rupture of mycobacterial cells with zirconium beads in a Biospec Products Mini-Beadbeater, and Shal, J.S. et al., J. Clin. Microbiol. 33 (2), 322-328 (1995) describe the lysis of heat-inactivated mycobacterial cells with the lysogenic agent guanidinium thiocyanate (GuSCN) and physical agitation with zirconium oxide beads. Also, U.S. Patent No. 5,374,522 describes methods of disrupting cells by applying ultrasonic energy to the cells in the presence of beads, and Hurley et al., Int. J. Systemic Bacteriology 38 (2):143-146 (1988) describe physical agitation of samples containing mycobacterial cells in the presence of distilled phenol and zirconium beads.

Alternatively, a solution containing mycobacterial microorganisms can be subjected to very intense ultrasonic bombardment in the presence of beads or particles which results in cell breakage. Typically, ultrasonic devices such as powerful ultrasonic probes (known as sonifiers or sonicators) are used in these processes. (See for example, Seiter, J. A. and Jay, J.M., "Application of Polyacrylamide Gel Electrophoresis to the Characterization and Identification of *Arthrobacter* Species," Int. J. Syst. Bacteriol., 30:460-465 (April, 1980)). However, significant amounts of heat are generated with high-powered probe devices of this type, and thus, cooling jackets or ice baths are required to reduce temperatures which can and often damage cellular nucleic acid. This damage to cellular nucleic acid from the high temperatures generated by use of sonicators with beads or particles has been shown by other researchers such as Salter, D.N. and Smith, R.H., "Protein Utilization in the Young Steer: Digestion and Nitrogen Retention of ¹⁵N-Labelled Rumen Bacterial Protein", British Journal of Nutrition, 51:531-539 (1984). These types of sonicating devices have measured outputs as high as 80-100 W as taught by Closs, O., et al., "The Antigens of *Mycobacterium bovis*, Strain BCG, Studied by Crossed Immunoelectrophoresis: A Reference System", Scand. J. Immunol., 12:249-263 (1980) and Allegro, H., U.S. Patent No. 3,558,066 entitled "Ultrasonic Extraction of Viable Antigens From Gram-Positive Bacteria," issued Jan. 26, 1971.

In contrast, U.S. Patent No. 5,374,522 teaches use of a sonicating device which operates at lower power densities. Thus, there is no need for cooling jackets or ice-baths

since the units lack the power to raise the temperature of the sonicating suspension to damaging levels. However, the method taught in this patent still requires the presence of beads or particles for effective cell lysis.

As recognized in U.S. Patent No. 5,374,522, rigorous physical grinding or shaking of organisms whether with or without beads presents considerable drawbacks. First, friction resulting from the physical interaction of grinding particles can create excessive heat which has deleterious effects on nucleic acids, and thus can render the nucleic acids unusable in subsequent hybridization procedures. Also, many of the organisms whose cells require such harsh conditions for extraction of cellular components are extremely pathogenic, and thus present health hazards when subjected to these physical manipulations in an open system. Also, the use of lysogenic chemical agents and/or enzymes such as SDS, GuSCN, proteinases, phenol/chloroform, etc. often adversely affects subsequent molecular biology processes for which the nucleic acids are accessed. For example, ionic and non-ionic detergents are known to inhibit nucleic acid amplification processes such as polymerase chain reaction (PCR) and strand displacement amplification (SDA), and carbon black which is commonly used to process glass beads is known to inhibit SDA.

SUMMARY OF THE INVENTION

The present invention provides an unexpectedly uncomplicated method for lysing cells to access nucleic acids therein. The present invention is a method for lysing cells wherein a sample containing cells is subjected to ultrasonic energy without the presence of beads. Following such lysis, nucleic acid from the cells is available for use in various molecular biology procedures.

BRIEF DESCRIPTION OF THE DRAWINGS

The various objects, advantages and novel features of the invention will be more readily appreciated from the following detailed description when read in conjunction with the appended figures in which:

- 5 Figure 1 is a graphical response surface plot showing the results of experiments to optimize the conditions of the method described herein.

DETAILED DESCRIPTION OF THE INVENTION

10 In a broad aspect, the method of the present invention provides a simple procedure for disrupting cells in order to access to cellular components including nucleic acids in a form useful for subsequent molecular biology procedures. The method's simplicity results from the limited manipulation of a sample and minimal additions to the sample.

15 The minimal additions to the sample results from the unexpected discovery that cells, including cells which are notoriously difficult to disrupt, such as mycobacterial cells, are disrupted by the application of ultrasonic energy in the absence of beads. Based on teachings such as those in U.S. Patent No. 5,374,522, it was believed by those skilled in the art that only with the inclusion of beads would the application of ultrasonic energy result in disruption of such difficult to lyse cells, such as mycobacterial cells.

20 A variety of ultrasonic baths are commercially available and useful for the present invention. Examples of suitable ultrasonic baths include those available from Branson Ultrasonics Corporation of Danbury, Connecticut, which markets a number of models under the Bransonic® name with tank capacities ranging from 0.5 to 5.5 gallons, and operating frequencies around 42 +/-6% kHz. Mettler Electronics® of Anaheim, Calif. also markets several models with tank capacities ranging from 2.1 quarts to 18 gallons. Blue Wave
25 Ultrasonics, Inc. of Davenport, Iowa markets a TT-0915 Table Top Cleaner with tank dimensions of 15.5"(l) x 9"(w) x 6"(d), and an operating frequency around a baseline 30 kHz. Lab-Line Instruments, Inc. of Melrose Park, Illinois offers 28 Lab-Line® ultrasonic bath

models with tank capacities ranging from 0.65 to 45 litres, including its model 9303 which operates at a fixed frequency of 35kHz, features industrial type transducers, and has a temperature range from ambient to about 80°C. Also, VWRbrand Aquasonic Ultrasonic Cleaners offer ultrasonic cleaners with bath volume capacities ranging from 2.5 to 7.5 gallons,
5 and operating frequencies of about 38.5 to about 40.5 kHz.

The Bransonic Ultrasonic Cleaner, Blue Wave Model TT-0915 and Aquasonic baths vary the ultrasonic frequency applied to the tank (*i.e.*, sweeping frequency). The resultant effect of sweeping frequency is the substantial reduction or elimination of standing waves and hot spots found when only a fixed frequency bath is used. The Bransonic Model 2510 is
10 particularly preferred for this sweeping frequency feature.

Such ultrasonic baths are recommended for cleaning tools, pens, jewelry, machinery, engine parts, nozzles, laboratory equipment, switches, locks, automobile parts, glass, ceramics, metals, hard plastics, etc. Ultrasonic cleaning baths such as these utilize a piezoelectric transducer such as, for example, lead zirconate titanate or barium titanate or a
15 magnetorestrictive transducer to convert electrical input energy into high frequency ultrasonic energy. This mechanical energy, or vibration, is then coupled into and transmitted through the liquid contained in the cleaning tank.

The term ultrasonic refers to frequencies just above the range of human hearing, hence about 20 kHz. Alternatively, ultrasonic energy can be delivered directly to the solution or
20 suspension of cells through, for example, a transducer. A solution or suspension of cells or microorganisms in purified or unpurified form can be placed in, for example, a vessel or well or a series of vessels or wells composed of a material capable of transmitting ultrasonic energy. The well is either attached to or is in proximity to a suitable transducer or other device capable of translating input energy into ultrasonic energy. The cells can be placed
25 directly into the well or series of wells which act as sample holders, or, alternatively the cells can be placed in containers and submerged in liquid contained within the well. The well can be capped with a suitable closure to prevent leakage or aerosolization.

While the method by which ultrasound disrupts cells has not been fully elucidated, it is postulated that ultrasonic waves traveling through a liquid consist of alternate compressions and rarefactions. If the amplitude of the wave is high enough, a phenomenon known as cavitation is produced. Cavitation is the making and breaking of microscopic bubbles. As these bubbles or cavities grow to what is known as resonant size, they collapse instantly and violently in one compression cycle, producing high local pressure changes or perhaps 20,000 atmospheres. This mechanical shock, which is felt at a distance of a few microns, is responsible for cellular disruption in the case of the high power density instruments. (Alliger, *H. Ultrasonic Disruption*, reprinted from *American Laboratory*, October 1975.)

The cells which are subjected to ultrasonic energy in the ultrasonic bath may be any cells which are to be disrupted. More specifically, cells which contain nucleic acid for use in a subsequent molecular biology application, such as bacterial, viral, fungal and other nucleic acid containing cells can be subjected to the method of the present invention.

The cells will be in a second liquid, that is a liquid other than the liquid in the ultrasonic bath. Such second liquid may be the sample in which the cells are found or a liquid to which a cell sample is added. The second liquid containing the cells is held in a vessel such as a closed tube, and placed in the first liquid in the ultrasonic bath.

The cells of the organism to be lysed can be in H₂O, but also can be in suitable buffers such as Tris-buffered saline (50mM Tris-HCl, pH8.0), phosphate-buffered saline (50mM sodium phosphate, 150mM NaCl, pH8.0), polymerase chain reaction buffer (10mM Tris-HCl, pH8.8, 50mM KCl, 1.5mM MgCl₂), React6 (buffer name React6 is registered by Bethesda Research Labs) (50mM Tris-HCl, pH7.1, 50mM NaCl, 50mM KCl, 6mM MgCl₂), sodium phosphate (pH 5.0 to 12.0), Trizma 9.0 (sigma; Trishydroxyaminomethylamine), and detergents such as 0.5% Tween 20 and 0.5% Nonidet P-40. The addition of detergents or any other agent which reduces surface tension of the liquid in which cells are suspended, or subjecting such liquid to any condition which reduces such surface tension, enhances the effects of the application of ultrasonic energy. A liquid of a

lesser surface tension will cavitate more readily than a liquid of higher surface tension. Optionally, the heated sample can be centrifuged, making available the supernatant and pellet for subsequent use.

Although the application of ultrasonic energy to cells in the absence of beads was found to be effective in the present method for disruption of cells, the method may be enhanced by alkalization of the second liquid in which the cells are located. More specifically, if the second liquid is rendered alkaline, the cell disruption method of the present invention is enhanced. Any suitable means for alkalization of the second liquid can be used in the present invention. Examples of suitable alkalization agents include KOH, NaOH, Ca(OH)_2 , Ba(OH)_2 , NH_3 and Na_2CO_3 .

A further enhancement of the method of the present invention is the heating of the first liquid in the ultrasonic bath. Maximum cavitation is known to occur in pure water at a temperature of approximately 71°C as taught at the website for CAE Blackstone, 9 North Main Street, Jamestown, New York 14701 (www.caeblackstone.com). Sufficient heat to enhance the disruption of cells with ultrasonic energy in the absence of beads is about 65°C to about 75°C .

Heating of a sample may be accomplished by any suitable method. The heat range for disrupting the cells of a particular organism is readily obtainable by titrating different temperatures for different amounts of time against release of desired cellular components from the cells of an organism. The heating will lyse the cells of the organism with subsequent release of intracellular components. One limitation on the heating is that the particular intracellular component of interest not be susceptible to destruction by the heat. Suitable heating means include water baths, microwave ovens, convection ovens, forced hot air ovens, and the like.

The heating time required for exposing intracellular components in the sample generally ranges from about two minutes to about forty-five minutes. The amount of heat and time of heat is readily found by sampling a portion of the cells of the organism to be lysed and

examining for signs of lysis (e.g., detection of intracellular components), depending on the source from which the intracellular component is to be obtained.

The method of the present invention is most effective when the second liquid is of an alkaline pH, the first liquid is heated to about 65°C to about 75°C and the cells are subjected
5 to ultrasonic energy ranging from about 35 kHz to about 45kHz.

In the Examples below, the Applicants compared the method of the present invention to traditional methods for disrupting cells.

The following examples illustrate specific embodiments of the invention described in this document. As would be apparent to skilled artisans, various changes and modifications
10 are possible and are contemplated within the scope of the invention described.

EXAMPLE 1

15 Comparison of Methods for Disrupting *Mycobacterium tuberculosis* Cells Using Sonication Without Beads vs. Physical Agitation with Beads

The purpose of this Example was to determine if release of amplifiable nucleic acid from *M. tuberculosis* ("M. tb.") cells using a method employing ultrasonic treatment without
20 beads is equivalent to release of amplifiable nucleic acid from M.tb cells using a method employing 5 m/sec physical agitation in the presence of zirconium beads. The comparison was based on amplification of a specific extracted nucleic acid target.

25 MATERIALS

SAMPLE PROCESSING REAGENTS

22 Sputa: clinical samples
M.tb Negative NaLc pool
BBL Mycoprep NaLc-NaOH
30 BBL Mycoprep PO₄ Buffer
2ml Labcraft sample processing tubes
Eppendorf micro centrifuge
Buffered GuSCN inhibitor removal wash
Sample Diluent
35 0.25N KOH
Tb Neutralization Solution
Zirconium beads encapsulated in onion skin glass

M.tb H37Rv -70°C frozen stock DOM 110196
Lab-line® Ultrasonic Bath model 9303
Savant CellPrep Agitator
BD Lysolyser

5

AMPLIFICATION AND ASSAY REAGENTS

Oligonucleotide Devices (ODs) - microtiter plates with SDA amplification primers, SDA bumper primers, SDA fluorescence detector probes, dUTP and buffers dried in each well.

- 10 Enzyme Devices (ENDs) - microtiter plates with restriction endonuclease (BsoBI), polymerase (Bst), dCsTP, dATP, dGTP and buffers dried in each well (ODs and ENDs are more completely described in co-pending U.S. Patent Application Serial No. 08/964,020 filed on November 4, 1997, the disclosure of which is specifically incorporated herein by reference).

15

PROCEDURE

- Twenty-two sputa samples were processed to NaLc pellets by the procedure recommended by the CDC. Briefly, 5-10 mls of sputa was mixed with an equal volume of mucolytic decontaminating NaLc/NaOH, using BBL Mycoprep. The solution was allowed to stand for 15-20 minutes with occasional mixing. The solution was then neutralized with BBL Phosphate buffer to bring the final volume to 50 mls. The solution was then centrifuged at 3,000xg for 20 minutes. After decanting the supernatant, 2 mls of phosphate buffer was used to resuspend the NaLc pellet.

- Each NaLc sample was divided into two 500 µl aliquots and placed into each of two 2 ml sample processing tubes. All 500 µl aliquots were then spiked with enough M.tb H37Rv to represent 1.25 CFU in a final amplification reaction. A NaLc pool, previously validated as not containing M.tb, was included as a negative control. Samples were then processed by either of two nucleic acid extraction procedures.

- In a first control procedure, samples were treated by a control condition involving physical agitation with zirconium beads. Briefly, 500 µl of NaLc sample was washed with 1ml Buffered GuSCN inhibitor removal wash. After centrifugation at 12,200xg for 3 min., the supernatant was decanted. Three sample diluent washes followed. After decanting the supernatant of the final wash, one 3 mm glass ball along with one capsule consisting of zirconium beads in onion skin glass was added to each tube. The tubes were then recapped

and rendered non-infectious by heating in a forced air oven to 105 °C for 30 minutes in a BD Lysolyser. M.tb lysis proceeded by subjecting tubes to 5 m/sec physical agitation in the presence of zirconium beads in a Savant CellPrep Agitator. More detailed descriptions of the removal of inhibitors and the process for agitation with beads may be found in co-pending

5 U.S. Patent Application Serial Nos. 08/774,476 filed on December 30, 1996, 08/614,230 filed on March 12, 1996, and 08/963,934 filed on November 4, 1997, the disclosures of which are specifically incorporated herein by reference.

In the experimental condition, M.tb lysis occurred by subjecting the processed sample to high temperature, alkaline conditions, and ultrasonic energy. The 500 µl sample was

10 washed, as in the control condition, with 1ml Buffered GuSCN inhibitor removal wash. After centrifugation at 12,200xg for 3 min., the supernatant was decanted. Three sample diluent washes followed. After decanting the supernatant of the final wash, 100 µl of 0.25N KOH was added. The recapped tubes were vortexed and then placed in a degassed Lab-line® Ultrasonic Bath model 9303. The bath temperature was 67 °C at the start of fixed frequency

15 35kHz treatment. Following 30 minutes sonication, the bath temperature was 75 °C. A 600 µl volume of Tb Neutralization Solution was then added and the recapped tubes were vortexed. The samples were then rendered non-infectious after being placed in a forced air oven to 105 °C for 30 minutes in a BD Lysolyser. Spiked samples were then amplified and assayed in duplicate. The negative control was run as a single replicate.

20 Aliquots (~150 µl) of a sample from the tubes treated by both M.tb lysis procedures, were dispensed into each well of the ODs. The wells of the ODs were covered, and the ODs retained at room temperature for 20 minutes. The ODs were then uncovered, and incubated at 75°C for 10 minutes, while the ENDs were pre-warmed for 10 minutes to 52°C.

After the 10 minute incubation, 100 µl aliquots from each well of the ODs were

25 transferred (pipetted) to a corresponding well in the ENDs. The ENDs were then sealed with an adhesive cover, and introduced into a fluorescence reader instrument as described in co-pending U.S. Patent Application Serial No. 08/929,895 filed September 15, 1997, the disclosure of which is specifically incorporated herein by reference.

The fluorescence signal from the wells of the ENDs were monitored for 60 minutes.

30 MOTA units, derived from the integration of the fluorescent signal curve over time, were

used for determining positive and negative M.tb nucleic acid amplification. The results are shown below.

RESULTS

5

Sputum ID	M.tb Spiked 1.25 CFU/rxn	M.tb Lysis by Using Agitation with Beads MOTA	M.tb Lysis by Using Sonication Without Beads MOTA
551	+	2323 801	4487 4186
2879	+	53709 50886	309 248
2170	+	1966 2606	2298 2061
66	+	1440 232	27014 8279
2194	+	4120 249	19688 9164
2895	+	47602 14301	21728 3814
2075	+	4354 2412	35329 36942
83	+	37455 53158	30893 29242
2990	+	29129 8699	31034 37232
78	+	4434 1112	118 1260
877	+	17257 3940	31416 30390
7 Subpool	-	251	357
3140	+	13754 32867	46556 42993
3316	+	45048 39390	6223 6420
3593	+	1010 326	518 2293
3614	+	202 253	1369 1464
H36198	+	4574 4258	23657 16531
M9287	+	20170 2988	8836 7597
F45663	+	14379 20878	2993 977
SF 1300	+	1275 699	901 538
1338	+	21540 3902	2422 1535
1360	+	4679 5244	20601 25445
1389	+	3045 2212	276 320
7 Subpool	-	358	300
MEAN		13293	13355
S. D.		16852	14311

For each experimental set of 11 samples, the data was statistically analyzed through a randomized complete block design with subsampling model. ROC analysis was used to calculate the sensitivity rates for the two sample processing procedures. Lastly, distributions were compared by empirical comparison through the Kolmogorov-Smirnov test.

There were no significant differences seen between the means of the two sample processing procedures in either experimental set. From the ROC analysis, sensitivity curves appeared to be equivalent for the two procedures. The distributions of MOTA generated by the two procedures appear to be identical (p-value of first experimental set 0.6208 and 0.8603 of second experimental set).

CONCLUSION

No apparent differences in the effectiveness of *M.tb* lysis procedures exists between a control method employing agitation with zirconium beads and a novel method employing ultrasonic treatment without beads. The control method has been shown to provide greater than 90% sensitivity of acid fast bacilli (AFB) smear negative, culture positive *M. tb* clinical samples as evidenced by: (1) G. E. Pfyffer *et al.* Study for the direct detection of *M. tb* complex in respiratory specimens presented at ICAAC 1997 abstract D92 (91.7% AFB smear negative sensitivity); (2) G. Woods *et al.*, who analyzed respiratory specimens and found a resolved sensitivity of 100% for AFB smear negative *M. tb* clinical samples (American Society for Microbiology 98th General Meeting abstract C-301); and (3) an evaluation of AFB smear negative, culture positive *M. tb* freshly processed NaLc which resulted in sensitivity of 92% (T. Fort *et al.* Presented at Clinical Microbiology and Infection, 8th Meeting, vol.3, supplement 2, May 1997 abstract P677).

EXAMPLE 2

Optimization of Novel Cell Lysis Method Employing Sonication Without Beads

The purpose of this experiment was to optimize the conditions for the novel cell lysis method employing ultrasonic treatment without beads. Temperature was varied from ambient to the previously tested 67-75⁰C. Alkaline conditions were varied from neutral pH 7.25 to alkaline pH 13.6. Duration of ultrasonic treatment varied from no treatment to 30 minutes.

Face-centered experimental design was utilized to determine significant and optimal conditions.

MATERIALS

5

SAMPLE PROCESSING REAGENTS

BBL Mycoprep PO₄ Buffer

2ml Labcraft sample processing tubes

Eppendorf micro centrifuge

10 0.25N KOH, pH 13.6

0.000416N KOH, pH 10.46

deionized H₂O, pH 7.25

pH 7.99, Tb Neutralization Solution

pH 8.69, Tb Neutralization Solution

15 M.tb H37Rv -70°C frozen stock DOM 040798

Lab-line® Ultrasonic Bath model 9303

BD Lysolyser

AMPLIFICATION AND ASSAY REAGENTS

20 As listed in Example 1

PROCEDURE

Ninety 2 ml sample processing tubes were numbered 1 - 90. One ml of BBL Mycoprep PO₄ Buffer was aliquotted to each tube. Each tube was then spiked with a volume of M.tb H37Rv to give a final concentration of 1.5 CFU per amplification reaction. All tubes were then microfuged at 12,200xg for 3 min. and the supernatant was decanted. All tubes were rendered non-infectious in a BD Lysolyser for 30 min. at 105 °C.

Tubes numbered 66-90 were resuspended with 100 µl of the appropriate KOH solution as determined by the experimental design found in the results section. The tubes were then placed in a degassed, 67.5 °C prewarmed Lab-line® Ultrasonic Bath model 9303 and sonicated for the time determined by the experimental design.

Tubes numbered 26-65 were resuspended with 100 µl of the appropriate KOH solution as determined by the experimental design found in the results section. The tubes were then placed in a degassed, 41.9 °C prewarmed Lab-line® Ultrasonic Bath model 9303 and sonicated for the time determined by the experimental design.

Tubes numbered 1-25 were resuspended with 100 µl of the appropriate KOH solution as determined by the experimental design found in the results section. The tubes were then

placed in a degassed, 23.1 °C ambient Lab-line® Ultrasonic Bath model 9303 and sonicated for the time determined by the experimental design.

Immediately following ultrasonic treatment, 600 µl of Tb Neutralization Solution was added to each tube to bring the final pH of the amplification reaction to pH 8.69. For pH 13.6 / 250 mM KOH treated samples, pH 7.99 Tb Neutralization Solution was used. pH 8.69 Tb Neutralization Solution was used for all other tubes. After vortexing the samples, the reaction volumes were amplified and assayed for specific M.tb target nucleic acid as in Example 1. The results are shown below.

10 RESULTS

Tube Number	Target Temperature °C	pH / [KOH]	Duration (min.)	MOTA (mean of 5 replicates)
1-5	19-30	7.25 / none	0	6013
6-10	19-30	7.25 / none	30	11148
11-15	19-30	10.46 / 416 µM	15	32662
16-20	19-30	13.60 / 250 mM	0	20153
21-25	19-30	13.60 / 250 mM	30	55027
26-30	44-45	7.25 / none	15	27788
31-35	44-45	10.46 / 416 µM	0	12100
36-40	44-45	10.46 / 416 µM	15	39772
41-45	44-45	10.46 / 416 µM	15	15693
46-50	44-45	10.46 / 416 µM	15	41613
51-55	44-45	10.46 / 416 µM	15	24301
56-60	44-45	10.46 / 416 µM	30	41073
61-65	44-45	13.60 / 250 mM	15	54081
66-70	69-80	7.25 / none	0	14852
71-75	69-80	7.25 / none	30	12089
76-80	69-80	10.46 / 416 µM	15	18544
81-85	69-80	13.60 / 250 mM	0	30508
86-90	69-80	13.60 / 250 mM	30	53863

Statistical analysis to determine significant factors found that a t-value of 4.763 (significance cutoff = 0.05) was found for pH with increasing pH optimal for M.tb lysis.

15 Duration of ultrasonic treatment was also significant with a t-value of 3.01. Increasing duration of ultrasonic treatment to 30 minutes was found to be optimal with a plateau effect within the 15-30 minute timeframe. Response surface plots generated from the experiment

presented in Figure 1 show the combinatorial effect of increasing pH and duration of sonication. The R-squared value of the experimental response was 0.838.

CONCLUSION

The optimal factors for the lysis of M.tb have been found to be extreme alkaline conditions, for 30 minutes, when 35 kHz of ultrasonic treatment is used. The use of alkali conditions may alter the conformation of proteins found in cell membranes, thus facilitating the extraction of nucleic acid, when used in conjunction with ultrasonic energy. Using a frozen stock of M.tb H37Rv, temperature was found to be insignificant.

EXAMPLE 3

Confirmation of Temperature Requirements for Alkali Sonication Without Beads

The purpose of this experiment was to confirm temperature requirements for the lysis of freshly prepared M.tb H37Rv cells. The previous examples have employed a frozen stock of M.tb H37Rv and have determined that heating of an ultrasonic bath is unnecessary for the effective lysis of M.tb.

MATERIALS

SAMPLE PROCESSING REAGENTS

BBL Mycoprep PO₄ Buffer
2ml Labcraft sample processing tubes
Eppendorf micro centrifuge
0.25N KOH, pH 13.6
pH 7.99, Tb Neutralization Solution
Freshly prepared M.tb H37Rv
Lab-line® Ultrasonic Bath model 9303
BD Lysolyser

AMPLIFICATION AND ASSAY REAGENTS

As listed in Example 1

PROCEDURE

BBL Phosphate Buffer was pipetted into 26 tubes, 1 ml volume for each tube. Eight tubes were spiked with an appropriate volume such that 1.07 particles of freshly prepared M.tb H37Rv were present in an amplification reaction. The addition of freshly prepared M.tb

was repeated in separate tubes at 0.71 and 0.38 particles/rxn. Two tubes were left unspiked with M.tb as a sample processing control.

All tubes were then centrifuged at 12,200xg for 3 min. and the supernatant was decanted. M.tb was rendered non-infectious by subjecting the tubes to 105 °C for 30 minutes in a BD Lysolyser.

One half of the tubes, at each concentration of M.tb spiked, was resuspended with 100µl of 0.25N KOH and vortexed. The tubes were then placed in a degassed, 63.0 °C prewarmed Lab-line® Ultrasonic Bath model 9303 and sonicated for 30 minutes.

The remaining tubes, at each concentration of M.tb spiked was resuspended with 100µl of 0.25N KOH and vortexed. The tubes were then placed in a degassed, 25.2 °C prewarmed Lab-line® Ultrasonic Bath model 9303 and sonicated for 30 minutes.

All tubes were neutralized immediately following ultrasonic treatment with 600 µl pH 7.99 Tb Neutralization solution and vortexed. The reaction volumes were amplified in triplicate and assayed for specific M.tb target nucleic acid as in Example 1. Freshly prepared M.tb H37Rv was plated for colony count quantification. The results are shown below.

RESULTS

M.tb 1.90 CFU/ml	Lysis Temp. 63.0-73.6 °C	Lysis Temp. 25.2-32.9 °C
	TB MOTA	TB MOTA
N=	12	12
Mean	33683	19996
Stdev	26465	23288
C.V.	72	117
M.tb 1.28 CFU/ml	Lysis Temp. 63.0-73.6 °C	Lysis Temp. 25.2-32.9 °C
	TB MOTA	TB MOTA
N=	12	12
Mean	33764	18171
Stdev	22794	13062
C.V.	68	72
M.tb 0.70 CFU/ml	Lysis Temp. 63.0-73.6 °C	Lysis Temp. 25.2-32.9 °C
	TB MOTA	TB MOTA

N=	12	12
Mean	19646	9707
Stdev	22131	11128
C.V.	113	115
M.tb 0 CFU/ml	Lysis Temp. 63.0-73.6 °C	Lysis Temp. 25.2-32.9 °C
	TB MOTA	TB MOTA
N=	3	3
Mean	993	844
Stdev	13515	9230
C.V.	1362	1094

Two-way statistical analysis of variance has determined that high temperature (63.0-73.6 °C) alkaline sonication provides moderately higher efficiency of M.tb lysis as compared to ambient temperature alkaline sonication. The p-value for this determination was 0.0980.

5

CONCLUSION

Maximum cavitation in pure water is known to occur at a temperature of approximately 71 °C as taught at the website for CAE Blackstone, 9 North Main Street, Jamestown, New York 14701 (www.caeblackstone.com). Although Example 2 did not find that temperature was a significant factor for alkaline sonication M.tb lysis, one can surmise that the effects of freezing on the stock of M.tb H37Rv masked any benefit that higher temperatures would provide. Using freshly prepared M.tb H37Rv in this example allows the benefits of temperature facilitated maximum cavitation to be realized. In order to detect M.tb which is present in low numbers, as in an AFB smear negative, culture positive clinical sample, the parameters of temperature, pH, and duration of ultrasonic treatment must be optimal.

15

While the invention has been described with some specificity, modifications apparent to those with ordinary skill in the art may be made without departing from the scope of the invention. Various features of the invention are set forth in the following claims.

20

WHAT IS CLAIMED IS:

1. A method for disrupting cells comprising:
 providing a sonic bath comprising a first liquid;
 placing into said liquid a vessel comprising cells in a second liquid; and
 5 subjecting said cells to ultrasonic energy from said sonic bath of sufficient power and
 duration to cause disruption of said cells in the absence of beads.
2. The method of claim 1 wherein said second liquid is at an alkaline pH.
3. The method of claim 2 wherein the temperature of said first liquid is about 65°C to
 about 75°C.
- 10 4. The method of claim 3 wherein the cells are mycobacterial cells.
5. The method of claim 1 wherein said second liquid is exposed to a substance or
 condition which reduces the surface tension of said second liquid.
6. The method of claim 5 wherein the temperature of said first liquid is about 65°C to
 about 75°C.
- 15 7. The method of claim 6 wherein the cells are mycobacterial cells.
8. A method for disrupting cells by applying ultrasonic energy to a sample of cells in a
 liquid, wherein said liquid is exposed to a substance or condition which reduces the surface
 tension of said liquid.
9. The method of claim 8 wherein said liquid is contained in a vessel and said vessel is in
 20 a sonic bath comprising a second liquid.
10. The method of claim 8 wherein beads are present in said liquid.
11. The method of claim 8 wherein said liquid is at an alkaline pH.
12. The method of claim 9 wherein the temperature of said second liquid is about 65°C to
 about 75°C.
- 25 13. The method of claim 12 wherein the cells are mycobacterial cells.

5

P-4278 application

Response surface plot for MOTAx1000 Holding Temp at 75




FIGURE 1

DECLARATION FOR PATENT APPLICATION

As a below named inventor, I hereby declare that:

My residence, post office address and citizenship are as stated below next to my name.

I believe I am the original, first and sole inventor (if only one name is listed below) or an original, first and joint inventor (if plural names are listed below) of the subject matter which is claimed and for which a patent is sought on the invention entitled CELL DISRUPTION METHOD USING SONICATION, the specification of which is attached hereto unless the following box is checked:

☐ was filed on as Application Serial No. and was amended on

I hereby state that I have reviewed and understand the contents of the above identified specification, including the claims, as amended by any amendment referred to above.

I acknowledge the duty to disclose information which is material to patentability as defined in Title 37, Code of Federal Regulations, §1.56.

I hereby claim foreign priority benefits under Title 35, United States Code, §119(a)-(d) of any foreign application(s) for patent or inventor's certificate listed below and have also identified below any foreign application for patent or inventor's certificate having a filing date before that of the application on which priority is claimed:

NUMBER	COUNTRY	DATE FILED (Day/Month/Year)	PRIORITY CLAIMED	
			YES	NO

I hereby claim the benefit under Title 35, United States Code, §119(e) of any United States provisional applications(s) listed below:

APPLICATION NUMBER	FILING DATE

I hereby claim the benefit under Title 35, United States Code, §120 of any United States application(s) listed below and, insofar as the subject matter of each of the claims of this application is not disclosed in the prior United States application in the manner provided by the first paragraph of Title 35, United States Code, §112, I acknowledge the duty to disclose information which is material to patentability as defined in Title 37, Code of Federal Regulations, §1.56 which became available between the filing date of the prior application and the national or PCT international filing date of this application:

APPLICATION SERIAL NO.	FILING DATE	STATUS (Patented, Pending, Abandoned)

I hereby appoint the following attorney(s) and/or agent(s) to prosecute this application and to transact all business in the Patent and Trademark Office connected therewith:




NAME	REG. NO.	NAME	REG. NO.	NAME	REG. NO.
Richard J. Rodrick	26,985	Richard E. Brown	30,934	Alan W. Fiedler	33,690
David W. Highet	30,265	Robert M. Hallenbeck	31,364	Bruce S. Weintraub	34,277
John L. Voellmicke	30,348	Donna R. Fugit	32,135	Susan A. Capello	34,560
Eric M. Lee	30,471	Royal N. Ronning	32,529	Vincent A. Castiglione	34,795
Allen W. Wark	30,503	Nanette S. Thomas	33,310	Arthur D. Dawson	35,113
				Jeremy Lack	35,813


Address all telephone calls to David W. Highet, Esq. at telephone number (201) 847-5317.

EXPRESS MAIL LABEL NO. EM368495356

Address all correspondence to Richard J. Rodrick, Esq., Becton Dickinson and Company, 1 Becton Drive, Franklin Lakes, New Jersey 07417-1880.

I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

FULL NAME OF SOLE OR FIRST INVENTOR	
Oscar J. Llorin	7/30/98
FIRST INVENTOR'S SIGNATURE	DATE
	
RESIDENCE	
228 Glenmore Avenue, Catonsville, Maryland 21228 USA	
CITIZENSHIP	
USA	
POST OFFICE ADDRESS	
228 Glenmore Avenue, Catonsville, Maryland 21228 USA	
FULL NAME OF SECOND JOINT INVENTOR, IF ANY	
Matthew P. Collis	
SECOND INVENTOR'S SIGNATURE	DATE
	7/30/98
RESIDENCE	
Box 90-38, RR #1, Seven Valleys, Pennsylvania 17360 USA	
CITIZENSHIP	
USA	
POST OFFICE ADDRESS	
Box 90-38, RR #1, Seven Valleys, Pennsylvania 17360 USA	
FULL NAME OF THIRD JOINT INVENTOR, IF ANY	
Michael C. Little	
THIRD INVENTOR'S SIGNATURE	DATE
	7-30-98
RESIDENCE	
3 Lochwell Court, Baltimore, Maryland 21234 USA	
CITIZENSHIP	
USA	
POST OFFICE ADDRESS	

FULL NAME OF FOURTH JOINT INVENTOR, IF ANY	
James M. Harris	
FOURTH INVENTOR'S SIGNATURE	DATE
	7/30/98
RESIDENCE	
10942 Millbank Row, Columbia, Maryland 21044 USA	
CITIZENSHIP	
USA	
POST OFFICE ADDRESS	
10942 Millbank Row, Columbia, Maryland 21044 USA	

P-4278 declaration